

# MicroSnap<sup>™</sup> Coliform Method Validation for Dairy Products

#### Introduction

MicroSnap<sup>™</sup> technology is based upon the biochemical detection of analytes produced by growing and actively respiring bacteria. Specifically, the MicroSnap<sup>™</sup> Coliform assay is based upon the detection of the indicator enzyme, beta-galactosidase. This enzyme is produced and expressed by the bacteria intracellularly. The overexpression of this enzyme is prompted by the addition of inducers that have been added to the enrichment/incubation media. The enzyme is then measured using a pro-luciferin substrate that is cleaved by the producing light as the bacteria grow.

The bacteria express the beta-galactosidase throughout their growth cycle; therefore, enzyme expression is proportional to bacteria CFU levels. This value has been determined to be optimal following a 6 hour incubation at 37°C. Incubating beyond this timepoint to an 8 hour timepoint is used to produce presence-absence RLU levels and is considered a non-linear response. This study evaluates results from dairy product testing.

## **Equipment, Supplies and Reagents**

- MicroSnap<sup>™</sup> Coliform Kit (containing Incubation and Detection Devices)
- Incubator set at 37°C ± 2°C
- Buffered Peptone Water (BPW)
- Tryptone soy broth (TSB)

- Tryptone soy agar (TSA)
- Bacterial cultures
- Pipettes and tips
- Cell spreaders
- Vortex mixer

#### **Procedure**

#### **Starter Cultures Testing**

#### Part A: Protocol for Estimation of Backgrounds

Starter Cultures were prepared following the directions provided by an industry leading dairy manufacturer. The prepared Starter Cultures were then diluted as shown in the Table 1 (3 replicates (reps) per dilution and timepoint).

For each dilution, 1 mL of culture was added to a MicroSnap<sup>™</sup> Coliform Incubation Device and incubated at 37°C. Background detection of beta-galactosidase-like activity was performed at the start of incubation (T=0 hrs), and again after 6 hrs, 8 hrs and 24 hrs.

Table 1: Starter Culture Background Testing Plan										
Starter Culture Dilution	0 hours	6 hours at 37°C	8 hours at 37°C	Overnight (18 - 24 hours) at 37°C						
Neat	3 reps	3 reps	3 reps	3 reps						
1:10	3 reps	3 reps	3 reps	3 reps						
1:100	3 reps	3 reps	3 reps	3 reps						
1:1000	3 reps	3 reps	3 reps	3 reps						

#### Part B: Protocol for Estimation of Detection of Coliforms in Starter Cultures

Starter Cultures were prepared following the directions provided by the customer. Specifically, four overnight cultures were produced of the following organisms:

- a. Escherichia coli
- b. Citrobacter freundii
- c. Klebsiella oxytoca
- d. Enterobacter cloacae

Following incubation (24 hours), the four overnight coliform cultures were mixed, and a dilution series was produced as shown below in Table 2.

	Table 2	. Coliform Inocula	tion Plan	
Coliform Dilution Series	Estimated CFU/mL	Add 100 μL of each dilution to 1 mL Starter	Estimated Coliforms per mL in Neat Starter Culture	Replicate Incubation devices to be incubated at 37°C
-1	1 x 10 <sup>8</sup>	Culture, add		
-2	1 x 10 <sup>7</sup>	Incubation		
-3	1 x 10 <sup>6</sup>	Devices		
-4	1 x 10 <sup>5</sup>		10,000	x1
-5	1 x 10 <sup>4</sup>		1,000	x1
-6	1 x 10 <sup>3</sup>		100	x3
-7	1 x 10 <sup>2</sup>		10	x3
-8	1 x 10 <sup>1</sup>		1	x3
-9	1 x 10°		0.1	x3
-10	1 x 10 <sup>-1</sup>		0.01	x3
Blank	0		0	x3
			Total	x20

## Results Starter Culture 1

#### Part A: Estimation of Background

Starter Culture Background Estimation - Starter Culture 1										
Dilution	Device	Time (hours)				Brilliance <sup>™</sup> Agar: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation			
		0	6	8	24	0h	24h			
Neat	Blank	0	0	0	0	NG	NG			
Neal	+ve control	0	0	0	0	NA	+			
1,10	Blank	0	3	0	0	NA	NG			
1.10	+ve control	0	136	10,610	11,951	NA	+			
1.100	Blank	0	0	0	0	NA	NG			
1.100	+ve control	0	146	12,927	20,000	NA	+			
Note:	+ve control (E	. coli ATCC®	8739) inoculu	m of 64 CFU						

NG = no growth; + = confirmed growth. Brilliance<sup>™</sup> is a trademark of Thermo Fisher Scientific

## **Observations and Summary**

- 1. There was no background signal from Starter Culture 1.
- Signal inhibition was noted when the product was spiked with 64 CFU at each of the dilution level.
  With Neat product, the signal inhibition was the most. The Neat Starter Culture may be very inhibitory to the growth of the added Coliform.
- 3. Diluting the product 1:10 and 1:100 helped to reduce the signal inhibition further.

#### **Next Step**

The 1:10 Dilution was chosen for the protocol in Part B - Estimation of Detection of Coliforms in Starter Cultures.

Part B: Estimation of Detection of Coliforms in Starter Cultures

		Coliform D	etection - 1:1	0 Suspension	of Starter Cult	ture 1	
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation		
		0	6	8	24	Inoculum	6h
	Rep 1	3	20	457	5,567		
6	Rep 2	0	27	726	5,613	79	Confirmed growth
	Rep 3	0	16	620	11,225		
	Rep 1	2	5	178	6,760		Confirmed growth
6.5	Rep 2	0	11	436	i 4,876 41	41	
Rep 3	Rep 3	3	4	252	6,940		
	Rep 1	0	0	26	6,654		
7	Rep 2	2	2	112	7,197	7	Confirmed growth
	Rep 3	0	0	8	7,085		
	Rep 1	0	0	0	273		
7.5	Rep 2	0	0	0	34	5	Confirmed growth
	Rep 3	0	3	113	6,425		
	Rep 1	-	0	4	84		
8	Rep 2	2	0	0	236	0	Confirmed growth
	Rep 3	0	0	0	64		
	Rep 1	4	0	0	0		
Blank PDT + RPW	Rep 2	3	3	0	2	0	-
	Rep 3	3	0	3	0		
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	rowth	·	

- 1. Growth is evident in the 1:10 product dilution of Starter Culture.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. Limit of detection is 70 CFU with fractional detection at 8 hrs.

#### Part A: Estimation of Backgrounds

	Starter Culture Estimation of Backgrounds - Starter Culture 2										
Dilution	Device		Time	Brilliance <sup>™</sup> Agar: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation						
		0	6	8	24	0h	24h				
Neat	Blank	4	0	0	0	NG	NG				
Neat	+ve control      3      0      0      0	NA	+								
1.10	Blank	0	0	0	0	NA	NG				
1:10	+ve control	0	373	9,378	20,000	NA	+				
1.100	Blank	0	0	0	0	NA	NG				
1:100	+ve control	0	144	11,203	20,000	NA	+				
Note:	+ve control (E	. <i>coli</i> ATCC® 8	739) inoculum	of 64 CFU							

#### **Observations and Summary**

- 1. There was no background signal from Starter Culture.
- Signal inhibition was noted when the product was spiked with 64 CFU at each of the dilution levels. With Neat product the signal inhibition was the highest. The Neat Starter Culture may be very inhibitory to the growth of the added Coliform.
- 3. Diluting the product 1:10 and 1:100 helped to reduce the signal inhibition further.

#### **Next Step**

The 1:10 Dilution was chosen for the protocol in Part B - Estimation of Detection of Coliforms in Starter Cultures

Part B: Estimation of Detection of Coliforms in Starter Cultures

		Coliform D	etection - 1:10	0 Suspension	of Starter Cult	ture 2	
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation		
		0	6	8	24	Inoculum	6h
	Rep 1	2	110	1,106	13,485		
6	Rep 2	4	143	2,604	12,401	79	Confirmed growth
	Rep 3	3	282	2,376	12,377		
	Rep 1	3	70	1,347	20,000		Confirmed growth
6.5	Rep 2	3	83	2,683	14,914	41	
Rep 3	Rep 3	0	109	2,599	20,000		
Rep	Rep 1	3	17	282	13,763		
7	Rep 2	0	16	166	15,771	7	Confirmed growth
	Rep 3	0	8	126	14,095		
	Rep 1	2	15	165	14,255		
7.5	Rep 2	2	4	106	7,041	5	Confirmed growth
	Rep 3	2	12	138	14,853		
	Rep 1	0	0	0	1,013		
8	Rep 2	0	0	0	2,563	0	Confirmed growth
	Rep 3	2	0	0	357		
	Rep 1	2	0	0	0		
Blank PDT + BPW/	Rep 2	2	0	0	0	0	-
	Rep 3	0	0	0	0		
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	rowth		

- 1. Growth is evident in 1:10 product dilution of Starter Culture.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. Limit of detection is 50 CFU with detection at 8 hrs and fractional detection starting at 6 hrs.

#### Part A: Estimation of Backgrounds

	Starter Culture	Estimation o	f Backgrounds	- Starter Cult	ure 3	
Dilution	Device		Time (	(hours)		MRS: CFU/mL
		0	6	8	24	0h
	Rep 1	2	0	0	0	NA
Neat	Rep 2	2	0	0	0	NA
	Rep 3	0	0	0	0	NA
	Rep 1	0	0	0	0	NA
1:10	Rep 2	0	2	0	0	NA
	Rep 3	0	0	2	0	NA
	Rep 1	0	0	0	0	NA
1:100	Rep 2	0	0	0	0	NA
	Rep 3	0	0	0	0	NA
	Rep 1	0	0	0	0	NA
1:1000	Rep 2	3	0	0	0	NA
	Rep 3	0	0	0	2	NA
Blank	BPW only	0	0	0	0	NA

- 1. There was no background signal from Starter Culture.
- 2. The product was diluted 1:10 for the spiking study.

Part B: Estimation of Detection of Coliforms in Starter Cultures

		Coliform D	etection - 1:1	0 Suspension o	of Starter Cult	ture 3	
Dilution	Device		Time (	hours)		TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation
		0	6	8	24	0h	6h
4	Rep 1	0	12,678	20,000	NA	14,750	Confirmed growth
	Rep 1	0	2,765	20,000	NA		
5	Rep 2	0	2,435	20,000	NA	1,475	Confirmed growth
	Rep 3	0	2,762	20,000	NA		
	Rep 1	0	76	12,833	20,000		
6	Rep 2	0	154	12,916	20,000	148	Confirmed growth
	Rep 3	0	113	13,500	20,000		
	Rep 1	0	8	2,780	20,000		
7	Rep 2	0	7	2,656	20,000	15	Confirmed growth
	Rep 3	0	3	1,074	20,000		
	Rep 1	0	0	33	9,886		
8	Rep 2	0	0	2	0	2	Confirmed growth
	Rep 3	0	0	2	844		
	Rep 1	0	0	2	0		
9	Rep 2	0	0	0	0	<1	-
	Rep 3	0	0	0	0		
	Rep 1	0	0	3	3		
Blank PDT + RPW	Rep 2	0	0	0	0	0	-
	Rep 3	0	0	0	0		
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	owth		

- 1. Growth is evident in 1:10 product dilution of the Starter Culture.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. Limit of detection is 20 CFU with fractional detection at 8 hrs.

#### Part A: Estimation of Backgrounds

Starter Culture Estimation of Backgrounds - Starter Culture 4										
Dilution	Device		Time (	Brilliance <sup>™</sup> Agar: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation					
		0	6	8	24	0h	24h			
Neat	Blank	0	0	0	0	NG	NG			
Neat	+ve control	0	0	0	0	NA	+			
1.10	Blank	2	2	2	0	NA	NG			
1.10	+ve control	2	610	13,786	NA	NA	+			
1.100	Blank	0	0	3	2	NA	NG			
1.100	+ve control	0	297	20,000	NA	NA	+			
1.1000	Blank	0	3	4	0	NA	NG			
1:1000	+ve control	4	376	20,000	NA	NA	+			
Note:	+ve control (E	. coli ATCC® 8	3739) inoculu	m of 84 CFU						

- 1. There was no background signal from Starter Culture.
- 2. Signal inhibition was noted when the Neat product was spiked with 84 CFU. The Neat Starter Culture may be very inhibitory to the growth of the added Coliform.
- 3. Diluting the product 1:10 and 1:100 helped to reduce the signal inhibition further.
- 4. 1:10 dilution was chosen as the next step for spiking study.

Part B: Estimation of Detection of Coliforms in Starter Cultures

		Coliform D	etection - 1:10	0 Suspension o	of Starter Cult	ure 4	
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation		
		0	6	8	24	Inoculum	6h
4	Rep 1	0	10,041	12,089	13,989	1,3550	Confirmed growth
	Rep 1	3	1,280	13,512	13,489		
5	Rep 2	3	1,435	13,650	12,896	1,355	Confirmed growth
	Rep 3	0	1,418	13,056	11,345		
	Rep 1	2	33	7,327	13,056		
6	Rep 2	2	70	6,843	13,904	121	Confirmed growth
	Rep 3	2	52	5,980	11,198		
	Rep 1	0	7	185	20,000		
7	Rep 2	2	3	18	5,064	15	Confirmed growth
	Rep 3	0	9	883	20,000		
	Rep 1	2	0	3	0		-
8	Rep 2	0	0	96	1,274	5	Confirmed growth
	Rep 3	0	0	6	5,758		Confirmed growth
	Rep 1	0	5	3	0		-
9	Rep 2	2	4	0	3,771	1	Confirmed growth
	Rep 3	2	0	4	0		-
	Rep 1	2	3	4	2		
Blank PDT + RPW	Rep 2	NA	0	0	2	NG	-
	Rep 3	NA	0	4	2		
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	owth		

- 1. Growth is evident in 1:10 product dilution.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. Limit of detection is 50 CFU with fractional detection at 8 hrs.

#### Part A: Estimation of Backgrounds

Starter Culture Estimation of Backgrounds - Starter Culture 5										
Dilution	Device		Time (hours)				Brilliance <sup>™</sup> Agar Confirmation			
		0	6	8	24	Oh	24h			
Neat	Blank	0	0	0	0	NG	NG			
Neat	+ve control	3	0	0	0	NA	+			
1.10	Blank	2	3	3	2	NA	NG			
1.10	+ve control	0	843	11,186	NA	NA	+			
1.100	Blank	0	2	2	0	NA	NG			
1.100	+ve control	2	448	20,000	NA	NA	+			
1.1000	Blank	2	0	0	2	NA	NG			
1.1000	+ve control	0	384	20,000	NA	NA	+			
Note:	+ve control (E	. coli ATCC® 8	8739) inoculu	m of 84 CFU						

- 1. There was no background signal from Starter Culture.
- 2. Signal inhibition was noted when the Neat product was spiked with 84 CFU. The Neat Starter Culture may be very inhibitory to the growth of the added Coliform.
- 3. Diluting the product 1:10 and 1:100 helped to reduce the signal inhibition further.
- 4. The 1:10 dilution was chosen as the next step for spiking study.

Part B: Estimation of Detection of Coliforms in Starter Cultures

		Coliform D	etection - 1:10	O Suspension of	of Starter Cult	ure 5	
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation		
		0	6	8	24	Inoculum	6h
4	Rep 1	3	4,897	13,870	12,287	13,550	Confirmed growth
	Rep 1	2	1,017	12,413	20,000		
5	Rep 2	0	896	8,536	10,683	1,355	Confirmed growth
	Rep 3	3	1,002	8,651	13,060		
	Rep 1	0	45	2,214	13,707		
6	Rep 2	2	26	3,893	12,875	121	Confirmed growth
	Rep 3	0	57	3,798	11,354		
	Rep 1	2	5	672	12,025	15	Confirmed growth
7	Rep 2	0	7	191	13,800		
	Rep 3	0	6	323	11,010		
	Rep 1	0	6	139	13,235		Confirmed growth
8	Rep 2	0	0	6	3	5	-
	Rep 3	0	5	2	0		-
	Rep 1	2	0	5	0		-
9	Rep 2	0	4	0	0	1	
	Rep 3	0	3	10	2		
	Rep 1	3	0	5	0		
Blank PDT + BPW	Rep 2	NA	0	0	3	NG	-
	Rep 3	NA	0	7	3		
Note:	Each inoc	ulated replicat	e was streaked	to confirm gr	owth		

- 1. Growth is evident in 1:10 product dilution.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. The limit of detection is 50 CFU with fractional detection at 8 hrs.

#### Part A: Estimation of Backgrounds

Starter Culture Estimation of Backgrounds - Starter Culture 6									
Dilution	Device		Time	(hours)	Brilliance <sup>™</sup> Agar: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation			
		0	6	8	24	0h	24h		
Neat	Blank	0	0	2	0	NA	NG		
Neat	+ve control	0	0	0	0	NA	NG		
1.10	Blank	2	0	3	0	NA	NG		
1.10	+ve control	3	65	5,323	10,507	NA	+		
1.100	Blank	NA	NA	NA	NA	NA	NA		
1:100	+ve control	NA	NA	NA	NA	NA	NA		
1.1000	Blank	NA	NA	NA	NA	NA	NA		
1:1000	+ve control	NA	NA	NA	NA	NA	NA		
Note:	+ve control ( <i>E. coli</i> ATCC <sup>®</sup> 8739) inoculum of 78 CFU								

- 1. There was no background signal from Starter Culture.
- 2. Signal inhibition was noted when the Neat product was spiked with 78 CFU. The Neat Starter Culture may be very inhibitory to the growth of the added coliforms.
- 3. Diluting the product 1:10 helped to reduce the signal inhibition further.
- 4. The 1:10 dilution was chosen as the next step for spiking study.

Part B: Estimation of Detection of Coliforms in Starter Cultures

Coliform Detection - 1:10 Suspension of Starter Culture 6									
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation				
		0	6	8	24	Inoculum	6h		
	Rep 1	0	11	1,815	13,379				
6	Rep 2	0	18	1,609	10,855	87	Confirmed growth		
	Rep 3	0	28	2,862	20,000				
	Rep 1	0	9	1,201	11,408				
6.5	Rep 2	3	7	904	20,000	57	Confirmed growth		
	Rep 3	0	5	1,188	20,000				
	Rep 1	0	2	119	11,101	9	Confirmed growth		
7	Rep 2	2	0	112	13,678				
	Rep 3	4	3	276	12,275				
	Rep 1	2	0	21	7,733		Confirmed growth		
7.5	Rep 2	0	0	5	2,781	1			
	Rep 3	4	0	5	3,090				
	Rep 1	0	0	5	1,831				
8	Rep 2	0	3	0	1,010	0	Confirmed growth		
	Rep 3	6	0	5	3,510				
	Rep 1	0	4	0	0				
Blank PDT + BPW	Rep 2	NA	0	4	0	NG	NG		
	Rep 3	NA	4	2	7				
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	rowth				

- 1. Growth is evident in 1:10 product dilution.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. The limit of detection is 10 CFU with fractional detection at 8 hrs.

#### Part A: Estimation of Backgrounds

Starter Culture Estimation of Backgrounds - Starter Culture 7									
Dilution	Device		Time	(hours)	Brilliance <sup>™</sup> Agar: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation			
		0	6	8	24	0h	24h		
Neat	Blank	0	0	0	0	NA	NG		
Neat	+ve control	0	0	2	0	NA	NG		
1.10	Blank	3	0	0	0	NA	NG		
1.10	+ve control	4	70	6,830	14,154	NA	+		
1.100	Blank	NA	NA	NA	NA	NA	NA		
1:100	+ve control	NA	NA	NA	NA	NA	NA		
1.1000	Blank	NA	NA	NA	NA	NA	NA		
1:1000	+ve control	NA	NA	NA	NA	NA	NA		
Note:	+ve control ( <i>E. coli</i> ATCC <sup>®</sup> 8739) inoculum of 78 CFU								

- 1. There was no background signal from Starter Culture.
- 2. Signal inhibition was noted when the Neat product was spiked with 78 CFU. The Neat Starter Culture may be very inhibitory to the growth of the added coliforms.
- 3. Diluting the product 1:10 helped to reduce the signal inhibition further.
- 4. The 1:10 dilution was chosen as the next step for spiking study.

Part B: Estimation of Detection of Coliforms in Starter Cultures

Coliform Detection - 1:10 Suspension of Starter Culture 7									
Dilution	Device		Time (	TSA: CFU/ 100 μL	Brilliance <sup>™</sup> Agar Confirmation				
		0	6	8	24	Inoculum	6h		
	Rep 1	5	14	2,584	20,000				
6	Rep 2	2	33	3,585	20,000	87	Confirmed growth		
	Rep 3	0	11	1,239	20,000				
	Rep 1	5	14	1,382	20,000				
6.5	Rep 2	2	16	1,069	20,000	57	Confirmed growth		
	Rep 3	0	11	908	20,000				
7	Rep 1	7	9	688	20,000	9	Confirmed growth		
	Rep 2	2	4	43	12,320				
	Rep 3	0	6	39	13,690				
	Rep 1	2	6	33	20,000		Confirmed growth		
7.5	Rep 2	4	5	3	6,109	1			
	Rep 3	0	8	34	12,156	_			
	Rep 1	5	4	0	15		-		
8	Rep 2	0	0	0	2,071	0	Confirmed growth		
	Rep 3	5	0	6	5,934	_	Confirmed growth		
	Rep 1	3	0	0	5				
Blank PDT + BPW	Rep 2	NA	0	0	0	NG	NG		
	Rep 3	NA	5	5	3				
Note:	Each inoc	ulated replicat	e was streaked	d to confirm gr	rowth				

- 1. Growth is evident in 1:10 product dilution.
- 2. The spiked dilutions are growing and being detected at 6 hrs and 8 hrs.
- 3. Unspiked product sample (Blank) has no signal at any of the time points.
- 4. Limit of detection is 10 CFU with fractional detection at 8 hrs.

## **Condensed Whey Testing**

In addition to Starter Cultures, the dairy manufacturing customer asked us to help validate MicroSnap<sup>™</sup> with other matrices, including whey and cream. What follows is information regarding those studies.

#### Procedure

1. Four overnight cultures were prepared of the following organisms:

a. Escherichia coli b. Citrobacter freundii c. Klebsiella oxytoca d. Enterobacter cloacae

- 2. After 24 hours (the next day), the coliform cultures were mixed together and a dilution series was prepared in PBS. Accurate CFU/mL counts from dilutions were determined by plating the dilutions on TSA.
- 3. A 1:5 dilution of condensed whey was made in BPW and plated on Brilliance<sup>™</sup> Agar and Violet Red Bile Glucose Agar (VRBGA).
- 1 mL of diluted (1:5) condensed whey product was added to the MicroSnap<sup>™</sup> Incubation Devices, followed by the addition of 100 µL of the coliform dilution series (as shown in the table below) and the devices were incubated at 37°C.
- 5. The detection assay was run at 0 hrs, 6 hrs, 8 hrs and 24 hrs as shown below.
- 6. From each enrichment/incubation, after 6 hours, a streak confirmation was done on Brilliance<sup>™</sup> agar/VRBGA.

1:5 Product Suspension of Condensed Whey LoD										
Dilution	Device	Time (hours)				TSA: CFU/mL	Brilliance <sup>™</sup> Agar Confirmation	VRBGA Confirmation		
		0	6	8	24	0h	61	า		
	Rep 1	0	44	3,941	2,999		+	+		
	Rep 2	0	42	3,334	3,862	556	+	+		
	Rep 3	0	41	3,505	2,449		+	+		
	Rep 1	0	4	276	1,595	55	+	+		
6	Rep 2	0	0	360	2,102		+	+		
	Rep 3	0	0	281	1,184		+	+		
	Rep 1	0	0	11	706		+	+		
7	Rep 2	0	0	25	527	10	+	+		
	Rep 3	0	0	20	329		+	+		
	Rep 1	0	0	0	0	1	-	-		
8	Rep 2	0	0	0	2		-	-		
	Rep 3	0	0	0	0		-	-		
Blank	Rep 1	0	0	0	0	0	-	-		

### **Limit of Detection Results**

- 1. Condensed whey product was diluted to 1:5.
- 2. Growth is evident in the 1:5 product dilution.
- 3. The spiked dilutions are growing and being detected at 8 hrs at 50 CFU.
- 4. The unspiked sample (Blank) is not being detected after 24 hours. This shows that this product does not have indigenous beta-galactosidase producing bacteria.

## Cream (Half & Half Cream) Testing

### Procedure

- 1. Product was tested Neat or diluted 1:5 and 1:10 in BPW.
- 2. A mixture of coliform culture was produced.
- 3. The product suspensions were spiked with the diluted coliform mix to attain CFU/mL levels as shown in the table below.
- 4. The spiked suspensions were then plated onto Brilliance<sup>™</sup> Agar.

## **Limit of Detection Results**

Half & Half Cream LoD									
Dilution	Device	Time (hours)			Post Dilution Inoculum (CFU/mL) Brillance <sup>™</sup>	TSA CFU/mL	Brilliance <sup>™</sup> CFU/mL		
		0	6	8	24	0h	0h	6h	
	Rep 1	0	0	108	552		2,300	TNTC*	
Neat	Rep 2	0	0	425	685	20		TNTC	
neat	Rep 3	0	0	180	422	20		TNTC	
	Blank	0	0	0	0			TNTC	
	Rep 1	0	3	1,508	13,133			TNTC	
1.5	Rep 2	0	5	1,820	10,615	10		TNTC	
1.5	Rep 3	0	4	1,652	3,510			TNTC	
	Blank	0	0	0	0			TNTC	
	Rep 1	15	13	3,139	20,000			TNTC	
	Rep 2	0	34	6,987	20,000	50		TNTC	
1.10	Rep 3	2	15	3,575	20,000	50		TNTC	
	Blank	0	0	0	0			TNTC	

\*TNTC = Too Numerous to Count

- 1. For cream, growth is evident in Neat 1:5 and 1:10 product dilutions.
- 2. The spiked dilutions are growing and being detected at 6 hrs to 8 hrs at 20, 50 and 500 CFU.
- 3. The unspiked samples (Blank) are not being detected after 24 hours. This shows that this product does not have indigenous beta-galactosidase producing bacteria.

#### **Summary**

#### **Starter Cultures**

RLU Growth of Positive Control Spike										
	Background Enzyme	Neat Product	1 in 10 Product (6 hours)	1 in 10 Product (8 hours)	Threshold RLU at Mean ± 6 std dev Blanks					
Starter Culture 1	NO	0 RLU	136 RLU @ 64 CFU	10,610 RLU @ 64 CFU	10.2					
Starter Culture 2	NO	0 RLU	373 RLU @ 64 CFU	9,378 RLU @ 64 CFU	9.6					
Starter Culture 3	NO	0 RLU	NA	NA	0					
Starter Culture 4	NO	0 RLU	610 RLU @ 84 CFU	13,786 RLU @ 84 CFU	8.2					
Starter Culture 5	NO	0 RLU	843 RLU @ 84 CFU	11,186 RLU @ 84 CFU	8.6					
Starter Culture 6	NO	0 RLU	65 RLU @ 78 CFU	5,373 RLU @ 78 CFU	13.1					
Starter Culture 7	NO	0 RLU	70 RLU @ 78 CFU	6,830 RLU @ 78 CFU	16.4					

Limit of Detection (CFU/mL)									
	6 hours	8 hours	24 hours	Threshold RLU at Mean ± 6 std dev Blanks					
Starter Culture 1	410 - 790 CFU	50 - 70 CFU (F)	<10 CFU	10.2					
Starter Culture 2	50 - 70 CFU (F)	50 - 70 CFU	<10 CFU	9.6					
Starter Culture 3	1,480 CFU	20 CFU (F)	20 CFU (F)	10*					
Starter Culture 4	1,210 CFU (F)	50 CFU (F)	<10 CFU (F)	8.2					
Starter Culture 5	1,210 CFU	50 CFU	50 CFU (F)	8.6					
Starter Culture 6	870 CFU	10 - 90 CFU (F)	<10	13.1					
Starter Culture 7	570 CFU	10 - 90 CFU (F)	<10 (F)	16.4					
All Starter Cultures	Log 2.9 (857 CFU)	Log 1.68 (48 CFU)	Log 1.23 (17 CFU)	11.0 RLUs					

Note:  $10^*$  is the assumed cutoff threshold RLU (Threshold RLU at mean  $\pm$  6 Std Dev for blanks) based on the background signal of the other Starter Cultures tested at time 0; (F) = fractional

#### **Results Demonstrate That:**

- The MicroSnap<sup>™</sup> system will detect Coliforms in the dairy Starter Cultures at 6, 8 and 24 hours.
- The limit of detection (LoD) average at 6 hours is 857 CFU (range of 50 to 1,480 CFU).
- The LoD average at 8 hours is 48 CFU (range of 10 to 90 CFU).
- The LoD average at 24 hours or less is 17 CFU (range of presence/absence to 50 CFU).
- The fractionality shows that the CFU added may go through attrition and be lower when growth and detection occurs; hence, the true LoD will be lower.
- The calculated thresholds across the Starter Cultures using mean blanks ± 6 std dev from Time 0. RLUs remains steady except Starter Culture 3 which had all 0 RLUs for blanks.

#### Whey and Cream Testing

Limit of Detection (CFU/mL)									
	6 hours 8 hours 24 hours								
Condensed Whey	2,780 CFU	50 CFU	10 CFU	10*					
Half & Half Cream	500 CFU	20 CFU	10 CFU	10*					

Note:  $10^*$  is the assumed cutoff threshold RLU (Threshold RLU at mean  $\pm$  6 Std Dev for blanks) based on the background signal of the other Starter Cultures tested at time 0; (F) = fractional

#### **Results Demonstrate That:**

- The MicroSnap<sup>™</sup> system will detect Coliforms in the condensed whey cultures at 6, 8 and 24 hours.
- The MicroSnap<sup>™</sup> system will detect Coliforms in the half & half cream cultures at 6, 8 and 24 hours when diluted at least 1:5 in BPW. In Neat samples, Coliforms were detected at 8 and 24 hours.
- For whey, the limit of detection (LoD) average at 6 hours is 2,780 CFU and at 8 hours is 50 CFU.
- For cream, the LoD average at 6 hours is 500 CFU and at 8 hours is 20 CFU.
- The LoD average at 24 hours (or less) is 10 CFU.
- The fractionality shows that the CFU added may go through attrition and be lower when growth and detection occurs; hence, the true LoD will be lower.
- The calculated thresholds across the whey cultures and cream cultures all had 0 RLUs for blanks.

## Conclusions

For all Starter Cultures, growth was evident in the 1:10 product dilution of Starter Culture. In addition, all spiked dilutions showed detectable growth at 6 hrs and 8 hrs of incubation, while unspiked product samples (Blanks) had no signal at any of the time points.

For condensed whey, growth was evident in the 1:5 product dilution in PBS. In addition, all serial dilutions to -7 showed detectable growth at 8 hours, while unspiked samples (Blanks) had no signal at any of the time points.

For half & half cream, growth was evident at 8 hours in Neat 1:5 and 1:10 dilutions. Furthermore, growth was detected at 6 hours in cultures diluted 1:10 in BPW.

The only variability observed for the seven Starter Cultures was in the limit of detection, but values were within 7-fold of each other (ranging from 10 to 70 CFU). In most cases, fractional detection of Starter Cultures was observed at 8 hours. For whey and cream, fraction detection was observed at 8 hours and in some cases, 6 hours.

These results demonstrate that MicroSnap<sup>™</sup> Coliform can be used as a simple, rapid, accurate method to measure bacterial levels. Therefore, when paired with EnSURE<sup>™</sup> Touch, MicroSnap<sup>™</sup> Coliform can be used as an alternative method to plating, providing bacterial counts in less than 8 hours, rather than 24 - 48 hours.